MGN0614

Valid from 25/04/2019 to 28/02/2022

Product Sheet (EN)

2

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4

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6

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8

Eigenschaften (DE)

10

PRODUCT SHEET			Expired since 28/02/2022			
IDENTIFICATION OF THE MEDICAL DEVICE	Type	Rechargeable battery				
	Commercial designation	Batterie médicale Molift Quick Raiser 14.4V 2Ah Fils				
	Reference	MGN0614				
	EAN	3660766470681				
	Brand	NX				
	Compatible / Original battery	Compatible				
	Packaging	Unitaire				
RECOMMENDED USAGE						
Follow the instructions and recommendations specific to each model, using the technical instructions and document resources from the devices in which the battery is used						
Brands	Equipment	Models	PN			
Molift	Léve malade	Quick raiser	A0541000 / 0541000			
Identification	GENERAL TECHNICAL CHARACTERISTICS	Chemistry	NiCd			
		Type	SC HD 12S1P			
		IEC designation	12 KRM 23/43			
		Rated voltage	14.4V			
		Nominal capacity	2 Ah			
		Internal resistance Ω	<168 mΩ			
	The voltage and the actual capacity in use can be affected by several factors, especially the temperature, the discharge current, the pack's history (ex:use, storage), etc					
ELECTRICAL CHARACTERISTICS	CHARGE	Maximum charging voltage	20.4V			
		Standard charging current (15h)	200mA			
		Fast charging current (2,5h)	1000mA			
	DISCHARGE	Range of operating voltage	10.8V at 16.8V			
		Min tension in discharge	10.8V			
		Max discharge current	7000mA			
		Lifespan 80% DOD (0,5 C)	500 Cycles			
	MAINTENANCE	Frequency of maintenance charges at 20°C	6 Months			
	CONTROL ELECTRONICS	Electrical protection	Yes			
		Low voltage power cut	No			
		High voltage power cut	No			
		Max power cut voltage	Yes			
These devices not only designed to protect the pack in case of an equipment failure. They must not be used to control the discharge. The protection circuits have a response time of a few milliseconds.						
MECHANICAL CHARACTERISTICS	Dimensions (+/- 2mm)	Length	173 mm			
		Width	63 mm			
		Depth	33 mm			
		Weight (+/- 5g)	640 Gr.			
	Mechanical protection	shrink sleeve				
	Wire length (+/- 10mm)	110 mm				
	Terminal	Wires				
CONDITIONS OF USE, STORAGE, AND TRANSPORT	CONDITIONS OF USE	Charging temperature	0 at 45°C			
		Discharge temperature	-30 at 60°C			
	CONDITIONS OF STORAGE	Storage temperature	-20 at 60°C			
		Level of humidity	65,00 %			
		Max storage time	2 Years			
	TRANSPORT	UN code	Class 8			
			ADR/RID classification	Class 8		
			IMDG classification	Class 8		
			IATA classification	Class 8		

INSTRUCTIONS	COMMISSIONING	• Check the batteries and the connectors: wires not damaged, battery not swollen, burnt smell, oxidation of the connectors, leak... • Respect the polarity • Do a full charge with the adequate charger before the first use
	CHARGE	• Use an adequate charger • The battery is warmer during the charge: during the first charge, check that the battery's temperature stays in the temperature operating ranges. • In case of an abnormal heating, stop the charge by unplugging the charger within the realms of possibility, remove the battery from the equipment de l'équipement, have the equipment, the battery and the charger checked by a technician.
	CASE OF NON-WATERPROOF BATTERIES	• It is normal to observe a release of gas during the charge and use. Do not smoke. Place in suitable premises. • Open batteries need regular maintenance carried out by a qualified technician.
	CASE OF LITHIUM ION BATTERIES	There is a fire hazard with lithium ion batteries in the following cases: overload, short circuit, charge and use outside the voltage and temperature ranges.
	WARNINGS	• Read the instructions of your device. • Only use in compatible devices. • Respect the load and storage conditions. • Do not use if the battery is damaged, do not burn, do not pierce, do not dismantle or modify. The protection circuits protect the battery and the equipment: do not deactivate them.

Marking	CE
---------	----

EXPLANATION OF SYMBOLS	LABELLING	
		Catalogue reference
		Lot number
		manufaturer's address
		Conform to the directive CE 93/42 DM class I
		To recycle in a suitable salvage and recycling structure
		Read the product sheet and the instruction manual

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