PPD 26

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Thank you for purchasing your **PPD 26 pH Meter Thermometer**.

For best results with your device:

- **Read** these operating instructions carefully,
- **Observe** les précautions d'emploi

MEANINGS OF THE SYMBOLS USED

	WARNING, risk of DANGER! The operator must refer to this user's manual whenever this danger symbol appears.
	The CE marking indicates conformity with European directives, in particular LVD and EMC.
	The rubbish bin with a line through it indicates that, in the European Union, the product must undergo selective disposal in compliance with Directive WEEE 2002/96/EC. This equipment must not be treated as household waste.

CONTENTS


1. TECHNICAL CHARACTERISTICS.....	17
2. DESCRIPTION	18
2.1 TOP.....	18
2.2 REAR PANEL	19
3. USE.....	20
3.1 CONNECTION AND STARTING UP	20
3.2 ADJUSTING MEASURING THE TEMPERATURE	20
3.3 ADJUSTING AND MEASURING THE pH.....	20
3.3.1 <i>CALIBRATION AND ADJUSTMENT OF THE pH SCALE:</i>	20
3.4 MEASUREMENT OF POTENTIALS:.....	21
3.5 RECORDER OUTPUT (IN pH ONLY).....	22
3.6 REPLACING BATTERY	23
3.7 USING AN EXTERNAL POWER SUPPLY (OPTIONAL).....	23
4. VARIOUS POSSIBLE PROBLEMS	24
5. SAFETY RULES	25
6. MAINTENANCE.....	26
6.1 CLEANING.....	26
6.2 METROLOGICAL CHECK.....	26
6.3 REPAIRS	26
7. TO ORDER.....	27
7.1 DELIVERY CONDITION.....	27

1. TECHNICAL CHARACTERISTICS


2000-point liquid crystal display	
Ranges	0.00 to 14.00 pH -1999 to +1999 mV 0 to 150°C
Accuracy	± 1%
Manual temperature compensation	0 to 60°C
Temperature of use	0-50°C
Maximum altitude	2000m
Humidity	5 to 80% at 35°C
Pollution class	Category 2
Category of use	2
COMPLIANCE WITH STANDARDS:	
Safety	With the power supply provided: EN 61010-1. Normal environmental conditions (indoor use)
CEM	EN 61326, class B instrument, intermittent operation.
Power supply	9V power supply unit
Consumption	0,5W
Dimensions	187 x 106 x 54 mm
Weight	340g
Battery low indication	«  » symbol

2. DESCRIPTION

2.1 TOP

2.2 REAR PANEL

3. USE

3.1 CONNECTION AND STARTING UP

Unpack the instrument, the 9V battery, the KCl solution, and the operating manual.

Open the cover of the battery compartment. Insert the 9V battery, paying attention to the polarity markings.

Close the cover.

Connect the combined pH electrode or another combined electrode to the "pH" BNC connector on the back of the instrument.

Connect the temperature probe to the 3.5mm jack (T°C)

3.2 ADJUSTING MEASURING THE TEMPERATURE

Switch set to T°C

Temperature probe

(platinum resistance, 100Ω at 0°C) connected to the 3.5mm jack.

Display of the temperature

(the manual temperature adjustment knob is inoperative)

No temperature probe.

Adjust the display to the temperature of the medium to be analyzed using the manual temperature knob.

3.3 ADJUSTING AND MEASURING THE pH

3.3.1 CALIBRATION AND ADJUSTMENT OF THE pH SCALE:

3.3.1.1 Checks/precautions concerning the electrodes:

Before each pH measurement operation, make sure that:

- The glass bulb of the electrode does not contain any air bubble.
- The liquid electrolyte reference electrode is correctly filled with a potassium chloride solution of the right concentration (KCl, 1 mol.L⁻¹).
- The plug of the filling orifice of the reference electrode is open.
- There must be enough liquid and the sintered part of the reference must be immersed.

3.3.1.2 Calibration:

Example: Use of pH 4.00 and pH 7.00.

- Rinse the electrode in distilled water.
- Immerse the electrode in the pH7.00 buffer solution (always first).
- When the signal stabilizes, use the "Offset" knob to adjust the value displayed to 7.00
- Rinse the electrode in distilled water.
- Immerse the electrode in a pH 4.00 buffer solution.
- When the signal stabilizes, use the "Slope" knob to adjust the reading to 4.00
- Rinse the electrode in distilled water.
- Immerse the electrode in the medium to be analyzed.

3.4 MEASUREMENT OF POTENTIALS:

Set the switch to mV.

Example: Use of the Michaelis solution

The oxidation-reduction potential E of a solution is equal to the potential read (E read) plus the potential of the reference electrode (E ref).

For a reference electrode (Ag, AgCl, 1 mol.L⁻¹ KCl), E ref = 240mV at 20°C (Cf. table below).

For a Michaelis reference solution, E = 406mV.

Therefore E read = 406 – 240 = 166mV.

- Rinse the combined electrode or the separate electrodes (for example: a platinum electrode and a reference electrode) in distilled water.
- Immerse the electrode in the Michaelis solution.
- Use the "Offset" knob to adjust the display to +166 if necessary.
- Rinse the electrode in distilled water.
- Immerse the electrode in the medium to be analyzed.

Remark: the slope and temperature adjustment potentiometers are inoperative.

TYPE OF ELECTRODE	POTENTIAL
Normal hydrogen electrode	0
Calomel electrode: Hg – Hg ₂ Cl ₂ – KCl 1 mol.L ⁻¹ Hg – Hg ₂ Cl ₂ – KCl 1 satur�	+ 284 + 252
Silver, silver chloride electrode: Ag – AgCl – KCl 1 mol.L ⁻¹ Ag – AgCl – KCl 3 mol.L ⁻¹ Ag – AgCl – KCl 3,5 mol.L ⁻¹ Ag – AgCl – KCl 1 satur�	+ 240 + 215 + 208 + 204
Mercurous sulphate electrode: Hg – Hg ₂ SO ₄ – K ₂ SO ₄ mol.L ⁻¹	+ 666

Remark: The potential of the mercurous sulphate electrode depends on the nature of the liquid in contact (variation of a few mV).

3.5 RECORDER OUTPUT (IN pH ONLY)

Polarity of the recorder connector.

Correspondence:

pH	Recorder output
0	0V
7	1V
14	2V

3.6 REPLACING BATTERY

When the  symbol appears on the display unit, it is essential to replace the battery.

To do this, open the cover on the back of the instrument, then replace the spent 9V battery with a new battery, paying attention to the polarity.
Close the cover.

3.7 USING AN EXTERNAL POWER SUPPLY (OPTIONAL)

Output characteristics: 9V  100mA
2.1mm jack (« + » in centre)

4. VARIOUS POSSIBLE PROBLEMS

PROBLEM	ADVICE
No display	Check the condition of the battery, the polarity. Clean the terminals of the battery with a clean, dry cloth.
The pH fails to stabilize in all solutions measured	Check that the pH electrode is correctly immersed (reference junction, glass bulb). Check the connection of this electrode.
The instrument indicates a value close to 7.00 whatever the pH	Check the condition of the pH electrode (glass bulb). Replace the pH electrode.
During calibration in the pH 7.00 buffer solution, impossible to adjust the display to 7.00 Offset potentiometer at stop.	Check the level of the reference electrolyte. Replace the electrolyte of the 1 mol.L ⁻¹ KCl reference electrode (if liquid electrolyte).
During calibration in the pH 4.00 or pH 9.00 buffer solution, impossible to adjust the display to 4.00 and 9.00, respectively. Slope potentiometer at stop.	Check the position of the temperature potentiometer. Clean the electrode by dipping it a few times in a dilute hydrochloric acid solution (0.1 mol.L ⁻¹ HCl).
No display of the automatically measured temperature	Check the connection of the temperature probe.

5. SAFETY RULES

- The interior of the instrument must always be kept clean and dry.
- If the instrument is used in a way not in conformity with the specifications, the protection provided by the instrument may be impaired.
- Disconnect the instrument before opening the housing.
- The PPD 26 pH-meter is designed for people familiar with good laboratory practice. If the PPD 26 is not used in accordance with these operating instructions, the protection provided by the equipment may be impaired.
- The operation of the instrument can be perturbed by radiated electric fields or conducted emissions.
- The instrument must not be used in a context requiring permanent operation without human supervision.
- This instrument produces, uses, and can emit RF energy, and if not installed and used in accordance with the operating instructions, it can interfere with radio communications.
- In a residential area, the user of this equipment will probably cause interference, in which case the user must, at his/her own expense, do everything necessary to remedy the interference.

6. MAINTENANCE

The instrument contains no parts that can be replaced other than by trained and accredited personnel. Any unauthorized repair or replacement of a part by an "equivalent" may gravely impair safety.

6.1 CLEANING

Use a cloth moistened with clean water or a neutral detergent to wipe the instrument, then wipe with a dry cloth.

Do not use the instrument again until it is completely dry.

6.2 METROLOGICAL CHECK

Like all measuring or testing devices, the instrument must be checked regularly.

This instrument should be checked at least once a year. For checking and calibration, contact one of our accredited metrology laboratories (information and contact details available on request), at our Chauvin Arnoux subsidiary or the branch in your country.

6.3 REPAIRS

For all repairs before or after expiry of warranty, please return the device to your distributor.

7. TO ORDER

7.1 DELIVERY CONDITION

- 1 pH Meter Thermometer PPD 26
- 1 user manual
- 2 pH buffer solutions
- 1 9V power supply unit

All packed in a cardboard box.

FRANCE

Chauvin Arnoux Group

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75876 PARIS Cedex 18

Tél : +33 1 44 85 44 85

Fax : +33 1 46 27 73 89

info@chauvin-arnoux.com

www.chauvin-arnoux.com

INTERNATIONAL

Chauvin Arnoux Group

Tél : +33 1 44 85 44 38

Fax : +33 1 46 27 95 69

Our international contacts

www.chauvin-arnoux.com/contacts