ENTRETIEN

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Measure up

GB

Thank you for purchasing a **conductivity measurement cell**.
 For the best results with your probe:
 ■ **read** this user manual carefully,
 ■ **comply** with the precautions for use.

This cell is used to measure the electrolytic conductivity of solutions.

CE The CE marking indicates compliance with the European directives, particularly LV and EMC.

 The crossed-out dustbin means that, in the European Union, the product is the subject of selective disposal as per directive DEEE 2002/96/EC: this equipment must not be treated as household waste.

USE

- Connect the cell to the measuring instrument.
- Remove the protective cap.
- Immerse the poles in the solution to be measured.
- Rinse the cell with demineralized water after each measurement.

Measure the cell's constant during calibration by means of a reference solution with standard conductivity whose conductivity variation as a function of temperature is known (through a standard solution correspondence table) rather than simply trusting the indication provided by the manufacturer.

Indeed, careless use of the cell could change its geometrical characteristics, thus modifying the value of the constant. The value provided by the manufacturer is just a rough indication.

The constant is linked to the conductivity by the following formula:
 $\sigma = G \cdot k$

where σ is the solution's conductivity in S/cm
 G is the conductance in S
 k is the constant of the cell in cm⁻¹.

Each time after cleaning, check the constant of the cell by single-point calibration of your probe.

CLEANING

Conductivity cells do not require much maintenance but the poles must always be clean. After each measurement, rinse the cell with demineralized water and wipe the body with a clean, lint-free cloth without touching the poles.

Regular cleaning, calibration and correct cell storage conditions ensure accurate measurements.

To remove:

- calcareous deposits, use an acid cleaning product (0.1 mol/L sulphuric acid).
 - greasy films, use ethanol or propanol.
 - bacteria or algae, use a chlorinated product (e.g. Javel water (bleach solution)).
- Then rinse the cell with demineralized water as described above.

Do not use any products which attack glass (such as hydrofluoric acid) or PVC (such as acetone).

Do not touch the platinum poles with your fingers.

MAINTENANCE

Store your conductivity probe in the following conditions:

- Short-term storage (< 1 day): store in a dry place or in demineralized water.
- Long-term storage (> 1 day): store in a dry place.

Check the condition of your cell by making sure that there is still a sufficient layer of platinum present.

With time and repeated use, the platinum black on the poles may be altered. If so, it is possible to replatinize your probe (please contact us).

STATE AT DELIVERY

The conductivity cell is delivered in a cardboard box with:

- a shim,
- a protective cap,
- a bilingual user manual,
- a multilingual instruction sheet.

For accessories and replacement parts, please visit our website:

www.chauvin-arnoux.com

LIFE SPAN

Conductivity cells are subject to wear.

Some specific measurement environments (high-concentration acidic or alkaline environments, organic solvents, etc.) or excessive temperatures may reduce the product's life span and damage your cell.

As a result, the life span of your cell depends entirely on the way it is handled and the conditions for its use.

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